

Web Bluetooth Medical Dashboard

from physical hardware to digital data

Github repository





Table of Contents




- | | | | |
|----------|---------------------|----------|-----------------|
| 1 | System Architecture | 5 | User Interface |
| 2 | Technical Stack | 6 | System Workflow |
| 3 | Core Features | 7 | Page Views |
| 4 | Data with LINQ | 8 | Future Roadmap |

Our medical team: Dedicated professionals




George Khayat
22232005



Elia Ghazal
22330018



William Ishak
22232003



Bassam Farhat
22330047

Problem Statement

Platform Fragmentation	Installation Friction
Native apps per OS (iOS/Android/Windows) → more dev work + inconsistent UX.	Users must find, download, install → lower adoption (especially occasional use).



Solution: Web Bluetooth API

01.

Cross-platform compatibility.

02.

Zero-installation access via
URL.

03.

Real-time physiological
data acquisition.

System Architecture

Presentation Layer

- HTML5/Razor views
- CSS3
- jQuery for dynamic UI updates.


Application Layer

- ASP.NET Core 6 MVC controller handling logic and routing.

Data Layer

- In-memory storage processed via LINQ queries and exported via XML/JSON.

Technology Stack



Core Feature: Web Bluetooth Integration

Workflow:

1. User Gesture (Click)
2. `navigator.bluetooth.requestDevice()`
3. Connect to GATT Server
4. Subscribe to Characteristics (e.g., `0x2A37`)

```
// A. User Gesture (Click Event)
$('#connectHeartRate').on('click', async function() {

    // B. Request Device
    const device = await navigator.bluetooth.requestDevice({
        filters: [{ services: ['heart_rate'] }] // Standard Heart Rate Service
    });

    // C. Connect to GATT Server
    const server = await device.gatt.connect();

    // D. Subscribe to Characteristics
    const service = await server.getPrimaryService('heart_rate');
    const characteristic = await service.getCharacteristic('heart_rate_measurement'); // UUID:
0x2A37
    await characteristic.startNotifications();
    characteristic.addEventListener('characteristicvaluechanged', handleNotifications);
});
```

Data Processing With LINQ

We utilize C# LINQ (Language Integrated Query) for powerful in-memory data manipulation, ensuring type safety and code readability.


```
// 1. Filtering by Device Type (Where)
public List<HealthReadingDto> GetReadingsByDeviceType(string deviceType)
{
    return _readings
        .Where(r => r.DeviceType.Equals(deviceType,
StringComparison.OrdinalIgnoreCase))
        .ToList();
}

// 2. Real-time Aggregation (GroupBy)
public Dictionary<string, List<HealthReadingDto>> GetReadingsGroupedByDevice()
{
    return _readings
        .GroupBy(r => r.DeviceId)
        .ToDictionary(g => g.Key, g => g.OrderByDescending(r => r.Timestamp).ToList());
}

// 3. Calculating Statistical Averages (Average)
public Dictionary<string, double> GetAveragesByType()
{
    return _readings
        .GroupBy(r => r.DeviceType)
        .Where(g => g.Any())
        .ToDictionary(g => g.Key, g => g.Average(r => r.Value));
}
```

User Interface & Experience





Real-Time Feedback

The dashboard employs jQuery animations to provide immediate visual feedback (pulse effects) when new data packets arrive.


Define the Pulse Animation (CSS)

```
@keyframes pulse {  
  0% { transform: scale(1); box-shadow: 0 5px 20px rgba(0,0,0,0.2); }  
  50% { transform: scale(1.05); box-shadow: 0 10px 40px rgba(0,0,0,0.3); }  
  100% { transform: scale(1); box-shadow: 0 5px 20px rgba(0,0,0,0.2); }  
  
.pulse-animation {  
  animation: pulse 0.5s ease-in-out;  
}
```

Trigger it with jQuery when data arrives (JS)

```
function handleHeartRateData(value) {  
  // Update text  
  $('#heartRateValue').text(value);  
  
  // Trigger Animation  
  $('#heartRateCard').addClass('pulse-animation');  
  
  // Remove class after animation finishes so it can be triggered again  
  setTimeout(() => $('#heartRateCard').removeClass('pulse-animation'),  
            500);
```





Responsive Design

Built with Bootstrap 5 and custom CSS gradients, the interface adapts seamlessly to various screen sizes, maintaining usability across desktop and mobile devices.

Bootstrap Grid System (HTML)


```
<div class="row mb-4">
  <div class="col-12 col-md-3">
    <div class="card shadow stat-card">    <!-- Content -->
      </div>
    </div>
  </div>
```

Custom Mobile Adjustments (CSS)

```
@media (max-width: 768px) {
  .stat-card .display-3 {
    font-size: 2.5rem; /* Smaller font for mobile */
  }

  .btn-group .btn {
    display: block;    /* Stack buttons vertically
    width: 100%;     */
    margin-bottom: 10px;
  }
}
```

How Does The System Work?



Main Dashboard View

Medical/Fitness Dashboard

Web Bluetooth Devices

Connect Heart Rate Monitor Connect Thermometer Connect ESP32 Device Disconnect All

No devices connected. Click a button above to connect via Web Bluetooth.

Heart Rate

Temperature

ESP32 Sensor

Statistics (LINQ Aggregations)

Total Readings: 0

Connected Devices: 0

Last Reading: --

Device Types: 0


Recent Readings (Last 20)

Timestamp	Device Type	Device ID	Value	Unit	Notes
-----------	-------------	-----------	-------	------	-------

Export Data


Export as JSON Export as XML View Analytics Clear All Data

© 2026 - Medical/Fitness Health Dashboard - Web Programming Final Project





Live Demo



Analytics Model & Covid Data Fetch View

Health Data Management

Total Records

30

Devices

1

Device Types

1

Last Update

22:32:36

Future Roadmap

Persistence

- Transition from in-memory storage to a persistent database (SQL Server) for long-term tracking.

Mobile Apps

- Develop companion native apps using .NET MAUI for broader hardware support.

Analytics

- Integrate Machine Learning for anomaly detection and health trend prediction.

Thank you very much!

